

Integrating LabVIEW and Python

Why LabVIEW?

- Seamless hardware integration
- Easy GUI development
- Rapid iteration & development
- Built-in parallelism
- Software engineering & productivity tools
- LabVIEW ecosystem
- Real-Time
- FPGA
- ...



Why Python?

- Thousands of mature, open source packages



NumPy

matplotlib

pandas




SQLAlchemy



Examples of using LabVIEW and Python together

Take advantage of thousands of Python packages



NumPy

matplotlib

pandas



TensorFlow



OpenCV



SQLAlchemy



Flask



Requests
http for humans



Take advantage of thousands of Python packages

- **NumPy**: the fundamental package for scientific computing with Python.
- **SciPy**: provides many user-friendly and efficient numerical routines.
- **scikit-learn**: machine learning in Python.
- **scikit-image**: image processing in Python.
- **matplotlib**: a Python 2D plotting library.
- **iPython / Jupyter**: interactive computing with Python.
- **pandas**: high-performance, easy-to-use data structures & data analysis tools
- **requests**: Python HTTP for humans.
- **SQLAlchemy**: database abstraction library.
- **Tensorflow**: machine learning platform

Integrate with anything that has a Python API




Microsoft Teams











Leverage additional data visualizations from Python


Fruit Counts by Year




My Donut Plot






Default violin plot



Rectangular box plot



Notched box plot



Reuse existing Python scripts

Co-workers



```
46 def create_testresults_stackedbargraph (pass_data, fail_data, error_data, target_image_filepath, labels, title):
47     ind = np.arange(len(pass_data)) # the x locations for the groups
48     width = 0.35 # the width of the bars: can also be len(x) sequence
49
50     plt.bar(ind, error_data, width, color='darkgrey')
51     plt.bar(ind, fail_data, width, bottom=error_data, color='lightcoral')
52     plt.bar(ind, pass_data, width, bottom=fail_data, color='mediumseagreen')
53     # You can find a list of colors here: https://python-graph-gallery.com/python-colors/
54
55     plt.ylabel('Number of Tests')
56     plt.xticks(ind, labels)
57     plt.title(title)
58     plt.legend(['Error', 'Fail', 'Pass'])
59     plt.savefig(target_image_filepath.replace('\\', '/')) # replace all instances of '\' with '/'
60
61
27 def create_donutplot (data, target_image_filepath, labels, title, colors):
28     # Labels is a string list. Data is a numeric list. Title is string. Colors is a string list of colors.
29     fig, ax = plt.subplots(figsize=(5, 5), subplot_kw=dict(aspect="equal"))
30     wedges, texts = ax.pie(data, wedgeprops=dict(width=0.5), startangle=-40, colors=colors)
31     bbox_props = dict(boxstyle="square,pad=0.3", fc="w", ec="k", lw=0.72)
32     kw = dict(xycoords='data', textcoords='data', arrowprops=dict(arrowstyle="-"),
33             bbox=bbox_props, zorder=0, va="center")
34
35     for i, p in enumerate(wedges):
36         ang = (p.theta2 - p.theta1)/2. + p.theta1
37         connectionstyle = "angle,angleA=0,angleB=()"#.format(ang)
38         kw["arrowprops"].update({"connectionstyle": connectionstyle})
39
40     ax.set_title(title)
41     ax.legend(labels)
42     plt.savefig(target_image_filepath.replace('\\', '/')) # replace all instances of '\' with '/'
43     # plt.show()
44     plt.close()
```

Internet

A screenshot of a web browser showing a Stack Overflow question. The browser tab is titled 'python - Creating a Boxplot with...'. The URL is 'https://stackoverflow.com/questions/44119653/creating-a-boxplot-with-matplotlib'. The page title is 'Creating a Boxplot with Matplotlib'. The question text says: 'I am using python 3 and jupyter notebook. I have a pandas dataframe that is structured like this:'. Below this is a code snippet for a pandas DataFrame with columns 'location' and 'price', containing data for Asheville on April 25. The user says: 'I am simply trying to create a boxplot for each location to show the range of prices among items in each location. When I ran the following code:'. Below this is another code snippet: 'import matplotlib.pyplot as plt', 'import numpy as np'. The page has a sidebar with 'Stack Overflow' logo, 'Home', 'PUBLIC', 'Tags', 'Users', 'Jobs', 'Teams', and a 'Learn More' button. At the bottom of the sidebar, it says 'By using our site, you agree to our terms of service'.



machine learning python code



Methods of calling Python from LabVIEW

Methods


- Python Node
- Python Integration Toolkit for LabVIEW by Enthought

Methods of calling Python from LabVIEW

Python Node


Python Node

- Supported in LabVIEW 2018 and later



Python Node

- Open session & specify Python version



Python Node

- Call Python module (.py) and function

C:\MyPythonScript.py

```
import math
```

```
def Add(a, b):  
    return a+b;
```

```
def ConcatenateStrings(str1, str2):  
    return str1 + str2;
```

```
def EuclideanDistance(point1, point2):  
    xDiff = point1[0] - point2[0];  
    yDiff = point1[1] - point2[1];  
    return math.sqrt(xDiff*xDiff + yDiff
```

2

module path C:\MyPythonScript.py

function name Add

Python Node

Close
Python
Session

error out

return data type 0

a DBL

b DBL

a + b


Python Node – Supported Data Types

LabVIEW Data Type	Python Data Type
Integers	<code>int</code>
SGL, DBL	<code>float</code>
String	<code>str</code>
Boolean	<code>bool</code>
Array	<code>list</code> or <code>NumPy array</code>
Cluster	<code>tuple</code>

Note LabVIEW 2019 adds support for NumPy arrays.

Python Node

- Close the session



Python Node – Getting Started


1. Install LabVIEW 2018 or later.
2. Install Python.
 - Install supported version (2.7, 3.6)
 - Install same bitness (32-bit, 64-bit) as LabVIEW
3. Locate your Python script and function
 - Install any Python packages required by the Python script
4. Write/Run LabVIEW VI that calls Python script/function using Python Node



Python Node – Getting Started

1. Install LabVIEW 2018 or later.
2. Install Python.
 - Install supported version (2.7, 3.6)
 - Install same bitness (32-bit, 64-bit) as LabVIEW
3. Locate your Python script and function
 - Install any Python packages required by the Python script
4. Write/Run LabVIEW VI that calls Python script/function using Python Node


<https://www.python.org/>



Python 3.6 <https://www.python.org/downloads/release/python-368/>
Python 2.7 <https://www.python.org/downloads/release/python-2716/>

Python Node – Getting Started

1. Install LabVIEW 2018 or later.
2. Install Python.
 - Install supported version (2.7, 3.6)
 - Install same bitness (32-bit, 64-bit) as LabVIEW
3. Locate your Python script and function
 - Install any Python packages required by the Python script
4. Write/Run LabVIEW VI that calls Python script/function using Python Node




If using Python 3.6, enable this checkbox in the installer

Python Node – Getting Started

1. Install LabVIEW 2018 or later.
2. Install Python.
 - Install supported version (2.7, 3.6)
 - Install same bitness (32-bit, 64-bit) as LabVIEW
3. Locate your Python script and function
 - Install any Python packages required by the Python script
4. Write/Run LabVIEW VI that calls Python script/function using Python Node

MyPythonScript.py

```
1 import sys
2 import cv2
3 import numpy as np
4 import os.path as op
5
```




Command Prompt


```
C:\Python36\Scripts>pip install numpy
```

Python Node – Getting Started





1. Install LabVIEW 2018 or later.
2. Install Python.
 - Install supported version (2.7, 3.6)
 - Install same bitness (32-bit, 64-bit) as LabVIEW
3. Locate your Python script and function
 - Install any Python packages required by the Python script
4. Write/Run LabVIEW VI that calls Python script/function using Python Node





Python Node Demos



Reuse Python Visualizations (matplotlib) Demo




Reuse Python API for Slack (Demo)



Python Node – Tips


- If Python function has a NumPy array parameter, then marshal your LabVIEW numeric array as a NumPy array (instead of a list)



LabVIEW 2019*

Python Node – Tips


- Debugging
 - Add token to <labview>\LabVIEW.ini to show the Python console while Python Node is running



A screenshot of a Notepad window titled 'LabVIEW.ini - Notepad'. The window displays the contents of the LabVIEW.ini file. The text is as follows:

```
File Edit Format View Help
IconEditor.MostRecentTab=00000002
IconEditor.Position="472,285,1448,756"
HierStateEdgeStyle=3
ShowConsoleForPythonNode=True
```

The line `ShowConsoleForPythonNode=True` is highlighted with a red rectangular box.




A screenshot of a command prompt window titled 'C:\Program Files (x86)\National Instruments\Shared\NIPythonInterface\.\niPythonHost.exe'. The window shows the output of the Python console, which is 'Saying hello...'.

Python Node – Tips

- Debugging
 - Add token to <labview>\LabVIEW.ini to show the Python console while Python Node is running

Tip Place a breakpoint on the Close Python Session to keep the console window open



Python Node – Tips

- Read/write to global variable by using a wrapper function

Global variable →

```
my_global_variable = "abc"
```

Wrapper function →

```
def read_my_global_variable():  
    global my_global_variable  
    return my_global_variable
```

Wrapper function →


```
def write_my_global_variable(value):  
    global my_global_variable  
    my_global_variable = value
```

Python Node – Tips

- Installing a set of packages

- Use `pip freeze > requirements.txt` to generate file that lists installed packages.

```
C:\Users\ahsu\AppData\Local\Programs\Python\Python36-32\Scripts>pip freeze > C:\Users\Public\Documents\requirements.txt
```



- Use `pip install -r requirements.txt` to install the list of installed packages


```
C:\Users\ahsu\AppData\Local\Programs\Python\Python36-32\Scripts>pip install -r C:\Users\Public\Documents\requirements.txt
```

Methods of calling Python from LabVIEW

**Python Integration Toolkit
(By Enthought)**


Python Integration Toolkit

- Supported in LabVIEW 2015 and later
- Price: \$749




Python Integration Toolkit

- Specify which Python installation to use



Python Integration Toolkit

- Open session & specify Python module (.py)



Python Integration Toolkit

- Call Python function


C:\MyPythonScript.py

```
import math
```

```
def Add(a, b):  
    return a+b;
```


```
def ConcatenateStrings(str1, str2):  
    return str1 + str2;
```

```
def EuclideanDistance(point1, point2):  
    xDiff = point1[0] - point2[0];  
    yDiff = point1[1] - point2[1];  
    return math.sqrt(xDiff*xDiff + yDiff
```



Python Integration Toolkit

- Close the session




Python Integration Toolkit – Getting Started

1. Install LabVIEW 2015 or later.
2. Install Python Integration Toolkit
 - Automatically installs Python and Canopy IDE
3. Locate your Python script and function
 - Install Python packages required by Python functions via Canopy IDE
4. Run LabVIEW VI that calls Python script/function using Python Integration Toolkit API



Python Integration Toolkit – Getting Started

1. Install LabVIEW 2015 or later.
2. Install Python Integration Toolkit
 - Automatically installs Python and Canopy IDE
3. Locate your Python script and function
 - Install Python packages required by Python functions via Canopy IDE
4. Run LabVIEW VI that calls Python script/function using Python Integration Toolkit API




Python Integration Toolkit – Getting Started

1. Install LabVIEW 2015 or later.
2. Install Python Integration Toolkit
 - Automatically installs Python and Canopy IDE
3. Locate your Python script and function
 - Install Python packages required by Python functions via Canopy IDE
4. Run LabVIEW VI that calls Python script/function using Python Integration Toolkit API


MyPythonScript.py

```
1 import sys
2 import cv2
3 import numpy as np
4 import os.path as op
5
```



Python Integration Toolkit – Getting Started

1. Install LabVIEW 2015 or later.
2. Install Python Integration Toolkit
 - Automatically installs Python and Canopy IDE
3. Locate your Python script and function
 - Install Python packages required by Python functions via Canopy IDE
4. Run LabVIEW VI that calls Python script/function using Python Integration Toolkit API




Python Integration Toolkit – Benefits

- Easy deployment of LabVIEW EXE and Python environment
- Includes Canopy Python IDE
- Native functions for the following:
 - Get/set Python global variables
 - Display Python console
 - Specify which Python installation to use
- Supports Python classes
- Variety of shipping examples



Easy deployment of LabVIEW EXE and Python environment

1. **[Toolkit] Export Python Runtime**
 - Tools»Python Integration Toolkit»Export Python Runtime
 - Choose Python packages to export
2. **[LabVIEW EXE Buildspec] Create EXE**
3. **[LabVIEW Installer Buildspec] Create installer**
 - Installs LabVIEW-built EXE
 - Installs Python Runtime to subfolder of installed application. Python Runtime is scoped just to the application.
 - Calls Python Runtime as post-install step




Includes Canopy Python IDE

- Analysis and development environment for scientists and engineers using Python
- Provides curated set of packages managed through a package manager GUI
- Tailored to needs and workflows of scientists, analysts, and engineers





Additional Functions

Get and set Python global variables using native functions





Display Python console for debugging



Specify path to Python.exe

Variety of Shipping Examples

- Face detection using OpenCV package
- Frequency burst detection
- Call a cloud service
- Publish data to the web with plot.ly
- Python-powered web dashboard
- Get and set Python global variables
- Capture Printed Output
- Communication between concurrent Python script and LabVIEW VI
- ...



Enthought Python Integration Toolkit Demos

- Face detection using OpenCV package
- Frequency burst detection
- Call a cloud service
- Publish data to the web with plot.ly
- Python-powered web dashboard
- Get and set Python global variables
- Capture Printed Output
- ...

Additional Use Cases

Concurrent communication between LabVIEW and long-running Python operation

How

- LabVIEW calls Python function that launches a long-running operation in a new thread... which can communicate data via global variable or Python queue.
- LabVIEW and the long-running Python function can communicate via Python queues and/or global variables

```
def launch_complicated_operation():  
    """ Called from LabVIEW to start the operation in a parallel thread """  
    global worker_thread  
  
    worker_thread = Thread(target=do_something_complicated)  
  
    # This line is important: the "daemon" flag instructs Python not to wait  
    # for the thread to finish before exiting. We set it to True, so that  
    # when we use Close Session from LabVIEW, Python will be sure to exit.  
    worker_thread.daemon = True  
  
    worker_thread.start()  
  
def do_something_complicated():  
    """ Simulates a long-running Python operation that periodically puts  
    information into the queue.  
    """  
  
    status_queue.put("Starting up...")  
  
    for idx in xrange(10):  
        # See if LabVIEW has requested that we stop.  
        if not keep_running:  
            status_queue.put("Exiting early...")  
            break  
  
        # Do work. In this case, we just sleep for 1 second.  
        status_queue.put("Working... {}".format(idx+1))  
        time.sleep(1)  
  
    status_queue.put("Done.")
```

Demo

```
def is_running():  
    """ True if the worker thread is running, False otherwise """  
    if worker_thread is not None and worker_thread.is_alive():  
        return True  
    return False  
  
def check_queue():  
    """ Get the latest item on the queue, or return an empty string """  
    try:  
        result = status_queue.get_nowait()  
        if result is None:  
            result = ""  
    except Empty:  
        result = ""  
    return result
```


Calling LabVIEW from Python

- https://github.com/ni/python_labview_automation

Launch LabVIEW

Set VI control values

Run VI

Get VI indicator values

Get error message

Close LabVIEW


```
from labview_automation import LabVIEW
lv = LabVIEW()
lv.start() # Launches the active LabVIEW with the listener VI
with lv.client() as c:
    control_values = {
        "DBL Control": 5.0,
        "String Control": "Hello World!",
        "Error In": {
            "status": False,
            "code": 0,
            "source": ""
        }
    }
    indicators = c.run_vi_synchronous(
        vi_path, control_values)
    print(indicators['Result'])
    error_message = c.describe_error(indicators['Error Out'])
lv.kill() # Stop LabVIEW
```

Additional Resources


Download Python demo code

- <https://github.com/allenh-ni/labview-python-demo>

Python data visualizations



Slack API



Additional Content

- Calling Python Class Methods Using LabVIEW Python Node

<https://knowledge.ni.com/KnowledgeArticleDetails?id=kA00Z0000019UFmSAM&l=en-US>

- Using Optional Arguments with Python Node in LabVIEW

<https://knowledge.ni.com/KnowledgeArticleDetails?id=kA00Z0000019WoASAU&l=en-US>

NIWeek Hands-On: Integrate LabVIEW with Python

- **When:** Wed, May 22nd, 10:30 – 11:30am
- **Where:** Meeting Room 18C
- **Speaker:** Danielle Jobe, VI Technologies

4G LTE 12:54 PM

Surveys

Title
Processing at the Edge: Why a Platform-Based Approach Is Ideal for the IIoT

Time
Tuesday, 1:00 PM - 2:00 PM

Speaker(s)
Nick Butler

Nick Butler

*1. Please rate the session content on the following

Overall Quality
- select one -

Technical Level
- select one -

Relevance to your job
- select one -

Relevance to published title and abstract
- select one -

Nick Butler

Navigation icons: Refresh, Back, Forward, Home, App Drawer

Before you go,
take the survey.

Stay Connected During and After NIWeek



ni.com/niweekcommunity



facebook.com/NationalInstruments



twitter.com/niglobal



youtube.com/nationalinstruments